



Certificate No.: **115319A** 

# **CERTIFICATE OF CONFORMITY**

# Test Pieces, 9.5mm Stainless Steel

Product Code	Product Detail
P8982-95	Metal Detectable Test Wand 9.5mm Stainless Steel

### Contaminant

Diameter	9.5mm			
Material	AISI Stainless Steel JIS SUS 316			
Grade	100			
Standard of Manufacture	ISO 3290 International			
Density	8.027 kg/ltr = 8.027 g/cm3			
Heat No.	46482			
Batch Number	812048			

Materials supplied from ISO 9001:2008 / Din Standard Quality Assured Suppliers.

## Chemical and Test Data of Contaminant

Carbon WT	С	<0.08%	0.032	Chrome %	Cr	16-18%	16.850
Sillicum	Si	<1.0%	0.290	Nickel %	Ni	10-14%	10.010
Manganese %	Mn	<2.0%	<2.0%	Molibdeno %	Мо	2-3%	2.050
Phosphorus %	Р	<0.045%	0.040	Hardness	HRc	25	25.000
Sulphur %	S	<0.03%	0.010	Hardness	HRc	39	39.00

UNCONTROLLED IF PRINTED PAGE(S) 1 OF 3











Certificate No.:

115319A

### Manufacture of Test Piece Declaration

We hereby declare that our Test Pieces are manufactured from either the following materials and confirm:

- Cast Acrylic Sheets are in conformity with requirements of EU regulation 1935/2004/EU on materials and articles intended to come in to contact with food. The used monomers and other raw materials meet the requirements of EU regulation No 10/2011
- PTFE Grade 400 Sheets supplied are tested in accordance to FDA 177.1550
- ACETAL COPOLYMER CC34031 GREY FDA / CG34096 BLUE FDA / CS33273 RED FDA / CL33402E
   YELLOW FDA / CJ34045 GREEN FDA
- The Colour concentrate meets the requirements for repeated Food Contact according to FDA 21 CFR 178.3297
- ACETAL RESIN Rods are in conformity with Directive 2002/72/EC and FDA Regulation 21 CFR 177.2470
- SEMI RIGID HOSING Are in conformity with Directive 1935/2004/EU and FDA Regulation 21 CFR 177.1500
- POLYPROPYLENE All monomers and additives are in conformity with Directive10/2011/EC and FDA 177.1520

# Test Piece Life Span

- The contaminant or sphere is secured in a Medium or Housing so has an indefinite life span. Life span of the contaminant or sphere all depends on the housing of the test piece
- Test Pieces life expectancy depends on a great extent of the use it receives. Furthermore, customers
  often replace Test Pieces long before they become worn out due to changes in mishandling and
  damage to the Test piece.
- If a Test piece becomes chipped, broken or unreadable we recommend new test pieces.

**UNCONTROLLED IF PRINTED** 









Certificate No.:

115319A

This document was prepared on behalf of Klipspringer Ltd and the information included is provided in good faith and to the best of our knowledge correct at the time of writing. Klipspringer Ltd offers the information within this document as a guide only, and it does not represent any guarantee of the prescribed products in the sense of legal guaranteed relations. Klipspringer Ltd accepts no liability resulting from this document. It is the responsibility of the end user to ensure the items purchased are suitable for the intended application.

The date on the Certificate of Conformity relates to when the contaminants or spheres were purchased. This cannot be altered as it allows traceability for Audit purposes.

Supplier	Klipspringer Ltd		
Address	Rynor House, Farthing Road, Ipswich, Suffolk, UK. IP1 5AP.		
Telephone	+44 (0) 1473 461 800		
Email	sales@klipspringer.com		
Website	www.klipspringer.com		

Sheena Britton Technical Compliance Manager Klipspringer Date of Issue 13-03-13

Authorised by S. Britton

Revision No. 001

Revised by

H. Smith

UNCONTROLLED IF PRINTED PAGE(S) 3 OF 3





